

Autosampler



Model: UC3265

1. Sample positions: 120 (cup)
2. Injection repeatability: full loop injection: $RSD_6 \leq 0.3\%$,
partial loop injection: $RSD_6 \leq 0.5\%$ (injection volume $\geq 3 \mu\text{L}$)
3. Sample residue: $\leq 0.01\%$ (specified needle washing procedure)
4. Operation control mode: local control and RS-232 control
5. Display: LCD, 4.3-inch true color, animated graphical interface
6. Control: 32-bit embedded microprocessor with ARM7 core