

Analytical Balance

SY Series 0.1/1mg Touch Screen Balance



- 5.7-inch colour touch screen and convenient graphic menu
- With powerful software SUPERTEX1.0, SY balance have functions such as weighing, parts counting, percent weighing, check weighing, animal weighing, totalization, density weighing, formula weighing, target weighing, etc.
- Multi-users management, can store 5000 sets of test data, can display the sensor temperature and automatic temperature compensation in real time to ensure the accuracy of the instrument.
- Excellent wireless sensing operation function, can set the operation of automatic correction and peeling function, to avoid contamination of the instrument by chemical reagent operation
- Electromagnetic balance sensor
- External/Internal calibration
- RS232 and USB output interface

Model	Max.Capacity	Readability	Pan Size	Dimension/N.W.	Shipping Size/G.W.	Calibration			
SY104	100g	0.1mg	90mm	230*393*350mm/8kg	520*380*500mm/10kg	External			
SY204	200g	0.1mg							
SY304	300g	0.1mg							
SY103	100g	1mg							
SY203	200g	1mg							
SY303	300g	1mg							
SY503	500g	1mg				Internal			
SY104C	100g	0.1mg							
SY204C	200g	0.1mg							
SY304C	300g	0.1mg							
SY103C	100g	1mg							
SY203C	200g	1mg							
SY303C	300g	1mg							
SY503C	500g	1mg							
Std Accessories	Calibration Weight (for external calibration balance), Adapter, Manual, Dust Cover								