



FA/JA Series 0.1/1mg Balance

- Large blue backlighting 120*55mm LCD screen
- With the sensor technology of German STL software, it supports multi-point calibration and has more stable performance
- The windshield of all around antistatic coating can effectively shield the electrostatic interference
- Electromagnetic balance sensor
- External/Internal calibration
- 4-level shockproof
- Counting function
- The built-in RS232 standard interface can connect to printer or computer
- Unit: carat, ounce and gram

Model	Capacity	Readability	Pan Size	Dimension/N.W.	Shipping Size/G.W.	Calibration
FA1004B	0-100g	0.1mg	Φ80mm	350*215*340mm/8kg	520*380*500mm/10kg	External
FA1204B	0-120g					
FA2004B	0-200g					
FA2104B	0-210g					
FA2204B	0-220g					
FA1004C	0-100g					Internal
FA1204C	0-120g					
FA2004C	0-200g					
FA2104C	0-210g					
FA2204C	0-220g					
JA1003N	0-100g	1mg	Φ80mm	350*215*340mm/8kg	520*380*500mm/10kg	External
JA1203N	0-120g					
JA2003N	0-200g					
JA2603N	0-260g					
JA3003N	0-300g					
JA5003N	0-500g					
JA1003C	0-100g					
JA2003C	0-200g					
Std Accessories	Calibration Weight (for external calibration balance), Adapter, Manual, Dust Cover					