

## YK Series 0.1/1mg Balance



- Large blue backlighting 120\*55mm LCD screen
- Powerful functions like basic weighing, counting, density weighing, check weighing, percentage and animal weighing, etc.
- With the sensor technology of German STL software, it supports multi-point calibration and has more stable performance
- The windshield of all around antistatic coating can effectively shield the electrostatic interference.
- Electromagnetic balance sensor
- External/Internal calibration

Model	Capacity	Readability	Pan Size	Dimension/N.W.	Shipping Size/G.W.	Calibration
YK1004B	0-100g	0.1mg	Φ80mm	360*2205*350mm/8kg	520*380*500mm/10kg	External
YK1104B	0-110g					
YK1204B	0-120g					
YK2004B	0-200g					
YK2104B	0-210g					
YK2204B	0-220g					
YK1004C	0-100g					Internal
YK1104C	0-110g					
YK1204C	0-120g					
YK2004C	0-200g					
YK2104C	0-210g					
YK2204C	0-220g					
YK1003N	0-100g	1mg	Φ80mm	360*2205*350mm/8kg	520*380*500mm/10kg	External
YK1203N	0-120g					
YK2003N	0-200g					
YK2603N	0-260g					
YK3003N	0-300g					
YK5003N	0-500g					
YK1003C	0-100g					Internal
YK1203C	0-120g					
YK2003C	0-200g					
YK2603C	0-260g					
YK3003C	0-300g					
YK5003C	0-500g					
Std Accessories	Calibration Weight (for external calibration balance), Adapter, Manual, Dust cover					