

YP series Electronic Balance



- The new design promotes the reaction speed and fixed-point speed.
- The unique metal base and mechanical construction ensure a high stability and increase the anti-interference performance.
- It has various functions such as: multi-units' conversion, parts counting, animal weighing, etc.
- RS232C output interface.

Model	Capacity	Readability	Pan Size	Dimension/N.W.	Shipping Size/G.W.	Power
YP1002	0~100g	10mg	Φ130mm	280*200*60mm/1.8kg	320*280*150mm/2.5kg	communication & direct current
YP2002	0~200g	10mg				
YP3002	0~300g	10mg				
YP4002	0~400g	10mg				
YP5002	0~500g	10mg				
YP6002	0~600g	10mg				
YP10002	0~1000g	10mg				
YP20002	0~2000g	10mg				
YP30002	0~3000g	10mg				
YP2001	0~200g	100mg				
YP4001	0~400g	100mg				
YP6001	0~600g	100mg				
YP8001	0~800g	100mg				
YP10001	0~1000g	100mg				
YP20001	0~2000g	100mg	Φ165mm	320*220*100mm/2.5kg	360*300*185mm/3.3kg	
YP30001	0~3000g	100mg				
YP40001	0~4000g	100mg				
YP50001	0~5000g	100mg				
YP60001	0~6000g	100mg	255*190mm	320*220*100mm/2.5kg	360*300*185mm/3.3kg	
YP80001	0~8000g	100mg				
YP 100001	0-10000g	0.1g	255*190mm	320*220*100mm/2.5kg	360*300*185mm/3.3kg	
YP 200001	0-20000g	0.1g				
YP1000	0~1000g	1g	Φ130mm	280*200*60mm/1.8kg	320*280*150mm/2.5kg	
YP2000	0~2000g	1g				
YP3000	0~3000g	1g				
YP 5000	0-5000g	1g	Φ165mm	280*200*60mm/1.8kg	320*280*150mm/2.5kg	
YP 6000	0-6000g	1g				
YP 10000	0-10000g	1g	255*190mm	320*220*100mm/2.5kg	360*300*185mm/3.3kg	
YP 20000	0-20000g	1g				
Std Accessories	Calibration Weight (optional), Adapter, Manual, Dust cover					