

Mini Centrifuge CF0201001



Features

1. Fixed speed of 5400rpm
2. Low noise ≤ 50 dB
3. High safety by low operation voltage and door open braking system
4. Enclosed design reduces the risk of liquid splashing into the equipment
5. Excellent dynamic balance keeps stable running even if the rotor is unbalanced on two sides
6. Click-on design allows for very easy rotor change without tools
7. Equipped with 2 interchangeable rotors of 1.5/2mL x 6 and PCR 8 x 2

Specifications

Model	CF0201001
Max. Speed [rpm]	5400
Max. RCF [x g]	2000
Input	100-240VAC, 50/60 HZ
Output	24VDC, 2A
Run time	Continuous operation
Driving motor	DC motor
Rotor capacity	0.2/0.5/1.5/2.0ml x 6
	0.2ml x 16PCR strips
	0.2ml x 2 PCR 8 strips
Noise level	≤ 50
Dimension	6.9" x 5.5",
Weight	1.85kg