



Mini Dry Bath Incubator with Heating Lid



Features:

- Heating Lid
- LED display for easy reading
- Temperature is up to 100°C
- Simultaneous display of set/actual time and temperature
- 3D digital wheel adjustment of speed and temperature
- High temperature precision is up to $\pm 0.5^\circ\text{C}$
- Programmable operation can store up to 5 programs
- Low operation voltage (12v) can be used directly with car charger
- Various applications with different blocks

Specifications:

Model	TC0401006
Temperature Control Range	RT + 5°C~100°C
Timing Range	0~99h59min
Temp. Control Accuracy	$\pm 0.5^\circ\text{C}$
Display Accuracy	$\pm 0.1^\circ\text{C}$
Heating Time (40°C to 100°C Temp.)	$\leq 15\text{min}$
Dimension (mm)	164 x 116 x 98
Net Weight	0.9kg

ACCESSORIES:

