

## R100 Mini Magnetic Stirrer



- LED display screen, which displays speed, mode, timing, etc.
- With smart chip, the speed can be set and saved, 2-modes of rotation and timing function available.
- Five different speeds are available.
- Two rotation options of forward/reverse mode and one-way mode.
- ABS plastic outer shell and anti-slip foot pads ensure stirring stability.
- Equipped with strut and electrode clip, which is used as stirred electrode stand.
- Laboratory routine use, can be carried in the field and on the spot, as well as in the incubator.
- Waterproof design.

### Specifications:

Speed range	Five-speed
Maximum mixing capacity	1000ml
Face diameter	100mm
Maximum operating current	150mA
Maximum power consumption	1.5W
Size of the electrode clip mount	Φ16mm (2 holes), Φ13mm (1 hole)
Electrode holder's Landing distance	0~100mm
Power adapter	DC12V Power adapter
Dimensions and Net Weight	100 x 160 x 48mm /0.6kg
Shipping Size and Gross Weight	230 x 220 x 110mm /1kg