



## Laboratory 10 x 0.4L Magnetic Hotplate Stirrer



(RoHS) (REACH)

### Specification

Model	<b>MS310</b>
Stirring position	10
Work plate material	Stainless Steel with silicone
Work plate dimension [LxW,mm]	470x205
Driving mode	Electromagnetic drive
Max.Stirring quantity[L]	0.4x10
Heating temperature[°C]	150
Speed range[rpm]	0-1200
Speed increment[rpm]	10
Stirring position distance[mm]	93
Max.Magnetic bar[L,mm]	30
Speed display	LED
Time display	LED
Economic mode	Yes
Adjustable reverse rotation	Yes
Operation type	Continuous/timing operation
Time range	0-99h59min
Voltage [VAC]	100-120/200-240
Frequency [Hz]	50/60
Heating Power [W]	500
Power[W]	580
Permissible ambient temperature[C]	5-40
Permissible relative humidity	80%
Protection class	IP30
Dimension [LxWxH,mm]	630x210x70
Weight[kg]	10