

## KingKong Mini MI0102001D



### Features

1. OLED display provides easy reading of the speed and time.
2. Touch key is more user-friendly and durable
3. Special black tempered glass set-up plate is waterproof, shatterproof; very easy to clean and dissipate heat.
4. The most advanced magnetic coil technology allows for wear-free drive without moving parts
5. Automatic reverse rotation every 30 seconds provides excellent mixing performance.
6. Max stirring quantity [H<sub>2</sub>O] of 1000ml
7. Speed range of 15-1500rpm

### Specifications

Model	MI0102001D
Stirring Positions	1
Working plate material	Tempered glass
Speed range [rpm]	15 - 1500
Display	OLED
Max. Stirring quantity [H <sub>2</sub> O, mL]	1000
Motor rating input [W]	5
Motor rating output [W]	3
Max. Stirring bar [mm]	25
Set-up surface [mm]	Φ100
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity	80%
DIN EN60529 Protection class	IP 65
Voltage [VAC]	100 - 240
Frequency [Hz]	50/60
Input power	6W
Dimension [W x D x H, mm]	117 x 12 x 180
Weight [kg]	0.3