

Temper Glass LED Digital Magnetic Stirrer



CE RoHS

5 INCH LED DIGITAL MAGNETIC STIRRER

Features:

- Max stirring quantity (H₂O) is 5000mL
- Wide speed range from 50 to 1500rpm
- Touch key is more smart and user-friendly
- Wear-free running provides longer life time
- LED display offers easy reading of the speed and time
- Pioneering and patented design of electromagnetism-driven mechanism
- No risk of demagnetization of the permanent magnet as the traditional design
- Automatic reverse rotation every 30 seconds allows excellent mixing performance
- No concern about the mechanical wear and noise of the motor as the traditional design
- Tempered glass working plate is waterproof; strong acid and alkali resistant and very easy to clean



SPECIFICATIONS:

Model	MI0102011
Set-up plate dimensions	5 Inch
Set-up plate material	Tempered glass
Stirring Mechanism	Electromagnetism drive
Max. stirring quantity [H ₂ O]	5000mL (3000mL at 1500rpm)
Speed range	50-1500
Speed display	LED
Time display	LED
Time range	0-99H99Min
Permissible ambient temperature [°C]	5 - 40 °C (in the operation)
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP65
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Power [W]	25
Dimensions	172x247x30mm
Weight [kg]	1.627

